Primers/Amorces

› Tumeur :

Gène	Séquence 5' - 3'	Référence GenBank [®]	Longueur de l'amplicon (pb)
<i>Glut1</i> (<i>Slc2a1</i>)	Fwd : TGC CTT GGA TGT CCT ATC TG Rev : ACC AGG GCC TAC TTC AAA GA	NM_011400.3	71
<i>Hprt</i> (<i>Hprt</i>)	Fwd : CAG GCC AGA CTT TGT TGG AT Rev : TTG CGC TCA TCT TAG GCT TT	NM_013556.2	147
<i>Ppia</i> (<i>Ppia</i>)	Fwd : AGC ATA CAG GTC CTG GCA TC Rev : TTC ACC TTC CCA AAG ACC AC	NM_008907.1	127
<i>Tweak</i> (<i>Tnfsf12</i>)	Fwd : GCT ACG ACC GCC AGA TTG GG Rev : GCC AGC ACA CCG TTC ACC AG	NM_011614.3	130
<i>Vegfa</i> (<i>Vegfa</i>)	Fwd : CAC CCA CGA CAG AAG GAG AG Rev : TCT CAA TCG GAC GGC AGT AG	NM_001025250.3	87

Dr AMIROUCHE A.
(2020/2021)

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Variation

Welcome to NCBI

The National Center for Biotechnology Information advances science and health by providing access to biomedical and genomic information.

[About the NCBI](#) | [Mission](#) | [Organization](#) | [NCBI News & Blog](#)

Submit

Deposit data or manuscripts
into NCBI databases

Download

Transfer NCBI data to your
computer

Learn

Find help documents, attend a
class or watch a tutorial

Develop

Use NCBI APIs and code
libraries to build applications

Analyze

Identify an NCBI tool for your
data analysis task

Research

Explore NCBI research and
collaborative projects

Popular Resources

[PubMed](#)[Bookshelf](#)[PubMed Central](#)[BLAST](#)[Nucleotide](#)[Genome](#)[SNP](#)[Gene](#)[Protein](#)[PubChem](#)

NCBI News & Blog

November 13 NCBI Minute: Resources
for next-gen sequence analysis

06 Nov 2019

On Wednesday, November 13, 2019 at
12 PM NCBI staff will present a webinar

NIH Biomedical Data Science Codeathon in
Pittsburgh, Jan 8-10

30 Oct 2019

NCBI is pleased to announce a
Biomedical Data Science Codeathon in

- ✓ All Databases
- Assembly
- BioCollections
- BioProject
- BioSample
- BioSystems
- Books
- ClinVar
- Conserved Domains
- dbGaP
- dbVar
- Gene
- Genome
- GEO DataSets
- GEO Profiles
- GTR
- HomoloGene
- Identical Protein Groups
- MedGen
- MeSH
- NCBI Web Site
- NLM Catalog
- Nucleotide
- OMIM
- PMC
- PopSet
- Probe
- Protein
- Protein Clusters
- PubChem BioAssay
- PubChem Compound
- PubChem Substance
- PubMed
- SNP
- Sparcle
- SRA
- Structure

- NCBI Home**
- Resource List (A-Z)**
- All Resources
- Chemicals & Bioassays
- Data & Software
- DNA & RNA
- Domains & Structures
- Genes & Expression
- Genetics & Medicine
- Genomes & Maps
- Homology
- Literature
- Proteins
- Sequence Analysis
- Taxonomy
- Training & Tutorials
- Variation

to NCBI

Center for Biotechnology Information advances science and health by providing access to genomic information.

[NCBI](#) | [Mission](#) | [Organization](#) | [NCBI News & Blog](#)

Submit

Submit data or manuscripts to public databases

Download

Transfer NCBI data to your computer

Learn

Find help documents, attend a class or watch a tutorial

Develop

Use APIs and code to build applications

Analyze

Identify an NCBI tool for your data analysis task

Research

Explore NCBI research and collaborative projects

Search

Popular Resources

[PubMed](#)

[Bookshelf](#)

[PubMed Central](#)

[BLAST](#)

[Nucleotide](#)

[Genome](#)

[SNP](#)

[Gene](#)

[Protein](#)

[PubChem](#)

NCBI News & Blog

November 13 NCBI Minute: Resources for next-gen sequence analysis

06 Nov 2019

On Wednesday, November 13, 2019 at 12 PM NCBI staff will present a webinar

NIH Biomedical Data Science Codeathon in Pittsburgh, Jan 8-10

30 Oct 2019

NCBI is pleased to announce a Biomedical Data Science Codeathon in

GenBank release 234 is available

Gene

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Create RSS Save search Advanced

- Gene sources
- Genomic
- Categories
- Alternatively spliced
- Annotated genes
- Non-coding
- Protein-coding
- Sequence content
- CCDS
- Ensembl
- RefSeq
- RefSeqGene
- Status
- Current
- Clear all
- Show additional filters

Tabular 20 per page Sort by Relevance

Send to:

See [KHSRP \(KSRP\) KH-type splicing regulatory protein](#) in the Gene database
[ksrp](#) in [Homo sapiens](#) [Mus musculus](#) [Xenopus tropicalis](#) [All 5 Gene records](#)

Search results

Items: 1 to 20 of 60

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See also 3 discontinued or replaced items.

<< First < Prev Page 1 of 3 Next > Last >>

Name/Gene ID	Description	Location	Aliases	MIM
<input type="checkbox"/> KHSRP ID: 8570	KH-type splicing regulatory protein [<i>Homo sapiens</i> (human)]	Chromosome 19, NC_000019.10 (6473102..6424811, complement)	FBP2, FUBP2, KSRP, p75	603445
<input type="checkbox"/> Khsrp ID: 16549	KH-type splicing regulatory protein [<i>Mus musculus</i> (house mouse)]	Chromosome 17, NC_000083.6 (57021049..57031507, complement)	6330409F21Rik, Fbp2, Fubp2, Ksrp	
<input type="checkbox"/> khsrp ID: 100170586	KH-type splicing regulatory protein [<i>Xenopus tropicalis</i> (tropical clawed frog)]	Chromosome 3, NC_030679.1 (130039451..130057575, complement)	VgRBP71, fbp2, fubp2, ksrp	
<input type="checkbox"/> khsrp.S ID: 399189	KH-type splicing regulatory protein S homeolog [<i>Xenopus laevis</i> (African clawed frog)]		XELAEV_18000829mg, VgRBP71, fbp2, fubp2, khsrp-b, ksrp	
<input type="checkbox"/> khsrp.L	KH-type splicing	Chromosome 3L,	XELAEV_18019123mg,	

clear

Gene

Gene

Advanced

Full Report ▾

Send to: ▾

KHSRP KH-type splicing regulatory protein [*Homo sapiens* (human)]

Gene ID: 8570, updated on 3-Nov-2019

Summary

Official Symbol KHSRP provided by [HGNC](#)**Official Full Name** KH-type splicing regulatory protein provided by [HGNC](#)**Primary source** [HGNC:HGNC:6316](#)**See related** [Ensembl:ENSG00000088247](#) [MIM:603445](#)**Gene type** protein coding**RefSeq status** VALIDATED**Organism** [Homo sapiens](#)**Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Primates; Haplorrhini; Catarrhini; Hominidae; Homo**Also known as** p75; FBP2; KSRP; FUBP2**Summary** The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization (Min et al., 1997 [PubMed 9136930]; Gherzi et al., 2004 [PubMed 15175153]).[supplied by OMIM, Apr 2010]**Expression** Ubiquitous expression in testis (RPKM 37.2), endometrium (RPKM 25.1) and 25 other tissues [See more](#)**Orthologs** [mouse](#) [all](#)

Go to [reference sequence details](#)

Genomic Sequence: NC_000019.10 Chromosome 19 Reference GRCh38.p13 Primary Assembly

Go to nucleotide: [Graphics](#) **FASTA** [GenBank](#)

Nucleotide

Nucleotide

Advanced

GenBank

Send to:

Homo sapiens chromosome 19, GRCh38.p13 Primary Assembly

NCBI Reference Sequence: NC_000019.10

[FASTA](#) [Graphics](#)

LOCUS NC_000019 11710 bp DNA linear CON 09-SEP-2019

DEFINITION Homo sapiens chromosome 19, GRCh38.p13 Primary Assembly.

ACCESSION [NC_000019](#) REGION: complement(6413102..6424811)

VERSION NC_000019.10

DBLINK BioProject: [PRJNA168](#)Assembly: [GCF_000001405.39](#)

KEYWORDS RefSeq.

SOURCE Homo sapiens (human)

ORGANISM [Homo sapiens](#)Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi;
Mammalia; Eutheria; Euarchontoglires; Primates; Haplorrhini;
Catarrhini; Hominidae; Homo.

REFERENCE 1 (bases 1 to 11710)

CONSRM International Human Genome Sequencing Consortium

TITLE Finishing the euchromatic sequence of the human genome

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FASTA

Send to:

Homo sapiens chromosome 19, GRCh38.p13 Primary Assembly

NCBI Reference Sequence: NC_000019.10

[GenBank](#) [Graphics](#)

>NC_000019.10:c6424811-6413102 Homo sapiens chromosome 19, GRCh38.p13 Primary Assembly

```

AGAGTGTCTCCGCGCCGTGTGGAGCGAGGCCCTTGTTCGCCGCTTGAGCCGCCGCCGCCGCCCTCC
TCAGCTTCAGCCTCCGCGCCAGGCCCGGCCCGCCGCGCCATGTCGGACTACAGCACGGGAGGACCCCG
CCCGGGCCCGCCGCCCGCCCGGGCGGGGGGGGGGAGCCGGAGGCGCCGGGGAGGCCCTCCGCCGGCC
CGCCAGCGCGGGGACCGGGCGGGCGGTCCCGCGCGCCCGCCCGGGCGGGGGTCCGCCGGGG
CCCCCTCTCAGCCACCCGGGAGGGCGCCCGGGAATCCGCAAGGACGCTTTCGCCGACGCCGTGCAGCGG
GCCCGCCAGCTGAGGAGCGCGAGGGCTGAGGGCGCCCTCCCGCGCCGCCAGCCGGGGAGGGGGC
GGGTCACGTGCGCGCGCGCGCAGAGGTCCCGAGGGGGCGGGGCCGTCCGGGGTACGTGGGCGG
TAGCCGGGGGCTGGCCCTCCCGCTGCGCCCTCGAGTGGGATCCCAACGGGTGACCCCTGCAGGGAGTA
GGGGCGCACTGGGGTCAACCAGGAGCGCGCCCTGGGTGTGCAGGCGCCCGGAGCTCGGAGACGCACTTT
TGGGACCCCCCTTAGCCCCCGCCCTTAGGACATGCCTGGAGGCAGGGTCCACCTTCCCGGGCGGAG
GGCCTTTACTCTCCCTCCATAACAGGCAAAAGGCGAGTGTTCGCCCAACCCAGTCCCATTTCTGTCC
CCTCTCTAGCAGAAAGCGACATTTGCGGGCTGGAGGACAGCTCCTTCTTTCAGTGGGAACTGGAT
GGGAACTCTCTAGGGCTGAGCACTTGGCGAAATGGGAAAGACAGCATGGGACTGGCTCTTGGATGC
TTTGGTGGCTGAAAGCGTGGTCTGTCTCTTTGGGGAAAGAAGGGAGGGTCTCCCACCCCGTGTATTG
AGGCGTAATGACAAGGAAGTGTACATGCTGCTAGGACTGACCTTCCAGCTCCATTTGCAGCCTCCCCA
TGGCCTGGAGGGGTGGGGTCCGTAGAAATGGGAGGACTTTGCCCGCGACCGGCTCTGGGCTGTCCCAAT
GCCTCCGGGGCTAGAGACTGGCTCCTCCAAAGTGGGAGGTCCCGTTAAGTGAACGCTGTGTTCTGGCC
CTCTGTCTCTCTCTCACCCCTGGAGTTGGGTGGGTCTTAGCGTGCATTTGGTCTTATTGCCCATGGAAC
CGTCTGATTCTGGGACAGCTTGAAGGGAGGAAGCTGGTCTGTCAAGTTGGTGGGGACAGGCCCCCTGTGG
GCATCCCATGAGTTTGGCTCTAGGGCAGCGTAGAAAAATGGGAGACTGGTGGAAAGGGAGATTCTTAACA
TGGTTTGGGTTTCTAAACCCAGACTCTGATCCTAAGTTTGGTGGGTTGAGGTTTCTTAAAGTTGCTCT
TCACGTGTCTGAGATGTAACCTTGTGCCATAGAGTCTTGGTTTTTGTGTTGTTGTTGTTTTTTGGAA
ACAGTCTCAATATCACCCAGGCGGGGTGCAGTGGCAGCATCTCAGCACACCACAGCCTCCACCTCTCCA
GTTGAAGTGTCTCTCTGCTCAGGCTCCCGAGTAGCTGGGATTATAAGTGTGCACCCACACAGCCAGCT
    
```

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Primer-BLAST

A tool for finding specific primers

Finding primers specific to your PCR template (using Primer3 and BLAST).

PCR Template

[Reset page](#)[Save search parameters](#)[Retrieve recent results](#)[Publication](#)[Tips for finding specific primers](#)

Enter accession, gi, or FASTA sequence (A refseq record is preferred)  [Clear](#)

NC_000019.10

Range

Forward primer From To
Reverse primer [Clear](#)

Or, upload FASTA file

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Primer Parameters

Use my own forward primer
(5'->3' on plus strand)

[Clear](#)

Use my own reverse primer
(5'->3' on minus strand)

[Clear](#)

PCR product size

Min Max

of primers to return

Primer melting temperatures
(T_m)

Min Opt Max Max T_m difference 

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section 

Exon junction span

Exon junction match

Exon at 5' side Exon at 3' side

Minimal number of bases that must anneal to exons at the 5' or 3' side of the junction 

Intron inclusion

Primer pair must be separated by at least one intron on the corresponding genomic DNA 

Amorces PCR

› Tumeur :

Gène	Séquence 5' - 3'	Référence GenBank [®]	Longueur de l'amplicon (pb)
<i>Glut1</i> (<i>Slc2a1</i>)	Fwd : TGC CTT GGA TGT CCT ATC TG Rev : ACC AGG GCC TAC TTC AAA GA	NM_011400.3	71
<i>Hprt</i> (<i>Hprt</i>)	Fwd : CAG GCC AGA CTT TGT TGG AT Rev : TTG CGC TCA TCT TAG GCT TT	NM_013556.2	147
<i>Ppia</i> (<i>Ppia</i>)	Fwd : AGC ATA CAG GTC CTG GCA TC Rev : TTC ACC TTC CCA AAG ACC AC	NM_008907.1	127
<i>Tweak</i> (<i>Tnfsf12</i>)	Fwd : GCT ACG ACC GCC AGA TTG GG Rev : GCC AGC ACA CCG TTC ACC AG	NM_011614.3	130
<i>Vegfa</i> (<i>Vegfa</i>)	Fwd : CAC CCA CGA CAG AAG GAG AG Rev : TCT CAA TCG GAC GGC AGT AG	NM_001025250.3	87